



Single Patient Drills

– safe, reliable and efficient



The best choice for you and your patients

Using Single Patient Drills provide safety, reliability and efficiency both when it comes to the surgical procedure and the handling in the clinic.

Features

Always sharp

Disposable and sterile

Uniform drill tip design

▶ Benefits

▶ Reduces the risk of overheating the bone due to blunt drills

▶ Eliminates the risk of contamination
▶ Less preparation saves time

▶ Simplified handling



Try Single Patient
Drills for your next
implant case


Convenient Drill Kit

The convenient and efficient Drill Kit is used for all implant cases. All drills are available in short and long versions for optimized handling in the anterior and posterior regions.



Uniform drill tip design

On the Single Patient Drills the tip height created by the point of the drill is 0.9 mm, irrespective of diameter.



Ask your sales representative for further information.

Astra Tech BioManagement Complex™

A successful implant system cannot be determined by one single feature alone. Just as with nature, there must be several interdependent features working together. The following combination of key features is unique to Astra Tech:

- **OsseoSpeed™** – more bone more rapidly
- **MicroThread™** – biomechanical bone stimulation
- **Conical Seal Design™** – a strong and stable fit
- **Connective Contour™** – increased soft tissue contact zone and volume

